

**iBloodGlucose361 is a blood glucose
log App
Operational Guidelines**

kutech

KUTECH INFORMATION CO.,LTD.

www.kutech.com.tw

iBloodGlucose361 Function v3.8.8

The iBloodGlucose361 is also a blood glucose log application. Browsing the data from notebooks in a time-consuming and labor-extensive manner is no longer required. There are 3 numerical ranges to measure with an empty stomach before meal: desirable lower than 100, IFG between 100-125, and diabetes higher than 125. There are 3 numerical ranges to measure with an two hours after a meal: desirable lower than 140, IGT between 140-200, and diabetes higher than 200. Blood sugar data maintenance: addition, correction, and deletion of blood sugar data are available. In addition, you can also record any notes about how you are feeling or other important information. KUO YUE always persists in giving consideration to both technology innovation and quality.

The software can have inexhaustible innovation and creation; therefore, it's enchanting. The recent technology develops in order to meet the user's needs. The renovation of the consumer software is one good example.

Lite Features :

- System develop based on iOS7 or later.
- Easy to use.
- Great user-interface
- To record your blood glucose(AC/PC), breakfast, lunch, dinner, bedtime and important information.
- Data transmission with e-mail
- **Export a text file** of your blood pressure data if you need a delimited text file to load into another program.
- **Export your blood pressure information directly to Excel (CSV)** if you're not sure what info you need and you'd rather delete and rearrange columns in a spreadsheet program.
- To display date range on charts header
- The charts shows the accumulated total of counts
- Computing and display bar charts statistics
- MAR function: Medication Administration Record.
- The statistics of ac and pc occurrences, percentages is made available.
- Blood glucose data maintenance: addition, correction, and deletion of blood glucose data are available.

iBloodGlucose361 Function v3.8.8

- There are 3 numerical ranges to measure with an empty stomach before meal: desirable lower than 100, IFG between 100-125, and diabetes higher than 125.
- There are 3 numerical ranges to measure with an two hours after a meal: desirable lower than 140, IGT between 140-200, and diabetes higher than 200.
- Adding a standard deviation statistic of blood glucose.
- Browsing the blood glucose reading were arranged according to months

iBloodGlucose361 Function v3.8.8

iPad 上午11:17 94%

Home

Blood glucose (6/6)

Feb,2014

READINGS	Blood Glucose Category	AC	PC (2 HRS)
All	Desirable	< 100	< 140
	IFG	100 - 125	
	IGT		140 - 200
AC	Diabetes	> 125	> 200

	Date	Readings	AC/PC	Meal	Notes	Edit
Desirable	02-23 08:30	97	AC	breakfast	Demo	i
IFG	02-23 13:33	122	AC	lunch	Demo	i
Diabetes	02-23 15:34	212	PC	lunch		i
PC (2 HRS)	02-23 19:35	142	AC	dinner	Demo	i
	02-23 21:36	324	PC	dinner	Demo	i
	02-24 11:16	145	PC	breakfast		i

iBloodGlucose361 Function v3.8.8

iPad 上午11:17 94%

[Home](#) 3 Numerical Ranges Monthly Statistics [calculate](#)

breakfast - + Feb - + 2014 - +

AC

3 Numerical ranges to measure

			<u>Percentages</u>	<u>Occurrences</u>
Desirable	< 100		-	-
IFG	100 - 125		-	-
Diabetes	> 125		-	-

Total :

PC (2 HRS)

3 Numerical ranges to measure

			<u>Percentages</u>	<u>Occurrences</u>
Desirable	< 140		-	-
IGT	140 - 200		-	-
Diabetes	> 200		-	-

Total :

iBloodGlucose361 Function v3.8.8

iPad 上午11:18 94%

< Home Standard Deviation Analysis calculate

Feb 2014

AC

Numerical ranges to measure	Meal Breakfast	Lunch	Dinner	Bedtime
n	-	-	-	-
Average	-	-	-	-
Std Dev	-	-	-	-

PC (2 HRS)

Numerical ranges to measure	Meal Breakfast	Lunch	Dinner	Bedtime
n	-	-	-	-
Average	-	-	-	-
Std Dev	-	-	-	-

Conclusions

Evaluation of Item	AC	PC
n	-	-
Estimated average glucose (eAG)	-	-
Standard deviation - mg/dL	-	-
Maximum Value- mg/dL	-	-
Minimum Value - mg/dL	-	-

iBloodGlucose361 Function v3.8.8

iPad 上午11:18 94%

< Home **Body Mass Index Services** calculate

Unit Imperial Metric

Your Height feet inches

Your Weight pound

4 Numerical ranges to measure		Level	Health Risk
Underweight	< 18.5	-	
Normal weight	= 18.5 - 24.9	-	
Overweight	= 25 - 29.9	-	
Obesity	>= 30	-	

Adults

Ideal Weight

-

pound

iBloodGlucose361 Function v3.8.8

iPad 上午11:18 94%

< Home New Data Entry Done

Blood Glucose

AC/PC AC PC

Meal Breakfast Lunch Dinner Bedtime

Notes

Fri Feb 21	8	15	
Sat Feb 22	9	16	
Sun Feb 23	10	17	
Today	11	18	AM
Tue Feb 25	12	19	PM
Wed Feb 26	1	20	
Thu Feb 27	2	21	